

# Epikut<sup>s</sup>



 **S.I.N.**  
Implant System  
**Abbreviated**

# EpiKut<sup>s</sup>



DOWNLOAD THE S.I.N. APP  
AND SEE IN AUGMENTED REALITY

**PLACE THE CELLPHONE CAMERA OVER THE IMAGE**



# Epikut<sup>s</sup> PLUS

EPIKUT S PLUS was idealized for you who wants to redefine the concept of dental implants. With a cutting and compressive design, double inverted support screws, combined with the best Plus surface, which accelerates the osseointegration, EPIKUT S PLUS is the perfect combination for those seeking superior results and with high predictability.



# THE UNBEATABLE COMBINATION OF DESIGN AND SURFACE THAT MAKES AN IMPLANT EPIC




## ➤ Indicated for all bone types

The exclusive macro geometry that features progressive cutting screws design makes EPIKUT S PLUS the state of the art for cases of immediate loading, low density bone, and post-extraction alveolus cases. Extremely versatile, EPIKUT S PLUS also allows its use in other clinical situations as long as the indicated drilling clinical protocol is followed.



## ➤ Accelerated bone healing

The high hydrophilicity, generated by an ultra-thin and homogeneous layer of hydroxyapatite, expands the activity of the proteins involved in the osseointegration process.



## ➤ Exclusive Plus surface

Developed in the main universities of Sweden, the Plus surface accelerates the osseointegration and promotes a higher bone quality, proven by over 50 preclinical studies.



## ➤ An implant with diverse possibilities

Morse Taper 16° connection making your clinical day-to-day easier.



## ➤ Clinical practicality

A single surgical kit for the installation of the complete EPIKUT S and EPIKUT S PLUS line.

# EPIKUT<sup>S</sup>

We recreated the concept of epic with EPIKUT S.

With a cutting and compressive design, double inverted support screws,  
this line provides more clinical practicality, predictability and high primary  
stability for those who seek superior results.




# THE NEW DEFINITION OF EPIC



## › Hybrid macro geometry, cylindrical body and conic apex

With an exclusive macro geometry and design of cutting screws, EPIKUT S is the best choice for cases of immediate load, low density bone and post-extraction alveolus, and it can also be used for all other clinical situations, always following the clinical steps suggested in its drilling system.



## › Double inverted support screws

Ensure greater primary stability and insertion torque.

## › Ultra-screwable

Profile of double and cutting screws ensure greater insertion speed of the implant.



## › Apex

Ease of insertion and bicortical fixation, especially in cases of very narrow remaining bone.



## › Exclusive cervical microthreads

Greater bone contact area and improves the dissipation of occlusal forces.



## › Adaptation accuracy

With exclusive and high stress resistant prosthetic components.


## › Manufactured in Cold Worked grade IV titanium

Super light metal, very resistant to corrosion, wear and fracture.



## › More options of prosthetic components for Morse Taper

Internal angulation of the EPIKUT S Morse Taper available at 16°.



## › Treatment on the entire surface

Double acid etching on the entire surface.




- Indicated for all types of bones, mainly for low density bones, post-extraction alveolar and immediate and/or late loading.
- It can be used for all other clinical situations, as long as the clinical steps suggested in the drilling system are followed.
- The exclusive macro geometry guarantees precision and agility at the time of surgery.

#### **INDICATIONS FOR CLINICAL USE:**

- 3.5 mm - Central incisors and lateral incisors
- 3.8 mm - Central incisors, canines and premolars
- 4.0mm - Premolars and molars
- 4.5 mm - Premolars and molars
- 5.0 mm - Premolars and molars

- 1.5 mm infra-bone installation
- Initial drill speed: 1200 rpm
- Speed of the drills 2.7 to 4.8mm: 800 rpm.
- Insertion speed: 20 to 40 rpm
- Maximum torque: 80 N.cm
- Immediate loading\*: recommended torque from 45 to 80 N.cm
- Includes cover screw of 2.0mm

\* Relative contraindication in patients with systemic or local problems and at the professional's discretion.



# EPIKUT S DRILLING SEQUENCE

## FOR SOFT TYPE BONES

Drilling sequence  
used for bone type  
IV.



		1.200 RPM		800 RPM						
	∅ DIAM. (mm)	FL 20 (A)	FHE 27 (B)	FHE 30 (C)	FHE 33 (D)	FHI 36 (E)	FHI 38 (E+)	FHI 40 (F)	FHI 43 (G)	FHI 48 (H)
ILM35xx	3.5	●		●						
ILM38xx	3.8	●		●	●					
ILM40xx	4.0	●		●	●	●				
ILM45xx	4.5	●		●	●	●	●			
ILM50xx	5.0	●		●	●	●	●	●	●	●

## FOR MEDIUM TYPE BONES

Drilling sequence  
used for bone type  
II and III.



		1.200 RPM		800 RPM						
	∅ DIAM. (mm)	FL 20 (A)	FHE 27 (B)	FHE 30 (C)	FHE 33 (D)	FHI 36 (E)	FHI 38 (E+)	FHI 40 (F)	FHI 43 (G)	FHI 48 (H)
ILM35xx	3.5	●		●	●	●				
ILM38xx	3.8	●		●	●	●	●			
ILM40xx	4.0	●		●	●	●	●	●	●	
ILM45xx	4.5	●		●	●	●	●	●	●	●
ILM50xx	5.0	●		●	●	●	●	●	●	●

- USE OF DRILL WITH COUNTERSINK FUNCTION - DEPTH OF 5 MM


## FOR HARD TYPE BONES

Drilling sequence  
used for bone type I.



	1.200 RPM				800 RPM					
	Ø DIAM. (mm)	FL 20 (A)	FHE 27 (B)	FHE 30 (C)	FHE 33 (D)	FHI 36 (E)	FHI 38 (E+)	FHI 40 (F)	FHI 43 (G)	FHI 48 (H)
ILM35xx	3.5	•	•	•	•					
ILM38xx	3.8	•	•	•	•	•	•			
ILM40xx	4.0	•	•	•	•	•	•	•	•	
ILM45xx	4.5	•	•	•	•	•	•	•	•	•
ILM50xx	5.0	•	•	•	•	•	•	•	•	•

## Technical measures



# PROSTHETIC SEQUENCE 16°

DIRECT SEQUENCE OVER THE IMPLANT

Single restorations



IMPLANT			
CODE EPIKUT S	CODE EPIKUT S PLUS	DIAM. (mm)	LENGTH (mm)
ILM 3585	ILM 3585N	3,5	8,5
ILM 3510	ILM 3510N	3,5	10
ILM 3511	ILM 3511N	3,5	11,5
ILM 3513	ILM 3513N	3,5	13
ILM 3515	ILM 3515N	3,5	15
ILM 3885	ILM 3885N	3,8	8,5
ILM 3810	ILM 3810N	3,8	10
ILM 3811	ILM 3811N	3,8	11,5
ILM 3813	ILM 3813N	3,8	13
ILM 3815	ILM 3815N	3,8	15
ILM 4085	ILM 4085N	4	8,5
ILM 4010	ILM 4010N	4	10
ILM 4011	ILM 4011N	4	11,5
ILM 4013	ILM 4013N	4	13
ILM 4015	ILM 4015N	4	15
ILM 4585	ILM 4585N	4,5	8,5
ILM 4510	ILM 4510N	4,5	10
ILM 4511	ILM 4511N	4,5	11,5
ILM 4513	ILM 4513N	4,5	13
ILM 4515	ILM 4515N	4,5	15
ILM 5085	ILM 5085N	5	8,5
ILM 5010	ILM 5010N	5	10
ILM 5011	ILM 5011N	5	11,5
ILM 5013	ILM 5013N	5	13
ILM 5015	ILM 5015N	5	15

## TITANIUM HEALING CAP

CODE	DIAM. (mm)	LENGTH (mm)
CIM 3502C	3.5	2.0
CIM 3504C	3.5	4.0
CIM 3506C	3.5	6.0
CIM 4502C	4.5	2.0
CIM 4504C	4.5	4.0
CIM 4506C	4.5	6.0

## PEEK HEALING CAP

CODE	DIAM. PLAT. (mm)	PROFILE DIAM. (mm)	HEIGHT (mm)
CPCM 0504	N/A	5	4
CPCM 0804	N/A	8	4
CPCM 0508	N/A	5	8

## OPEN TRAY TRANSFER

CODE	DIAM. (mm)
TMAIM 35C	3.5
TMAIM 45C	4.5

## ANALOG

CODE  
ANMP 3800



## DIGITAL ANALOG

CODE  
ADCM

## SCANBODY

CODE	DIAM. (mm)
JBSWCM	3.5



20 N.cm

**FINAL PRE-MILL  
BLANK**

CODE	DIAM. (mm)
AFTCM 1601 - DS	10
AFTCM 1602 - DS	14
AFTCM 1601 - MK	11.5
AFTCM 1602 - MK	15.8



20 N.cm

**TEMPORARY TITANIUM  
CYLINDER**

CODE	DIAM. (mm)	LENGTH (mm)
CPTM 3501 - H	3.5	1.0
CPTM 3502 - H	3.5	2.0
CPTM 3503 - H	3.5	3.0
CPTM 3504 - H	3.5	4.0
CPTM 4501 - H	4.5	1.0
CPTM 4502 - H	4.5	2.0
CPTM 4503 - H	4.5	3.0
CPTM 4504 - H	4.5	4.0



20 N.cm

**17° ANGLED CEMENTED  
ABUTMENT**

CODE	DIAM. (mm)	LENGTH (mm)
AIAM 3501C-H	3.5	1.0
AIAM 3502C-H	3.5	2.0
AIAM 3503C-H	3.5	3.0
AIAM 3504C-H	3.5	4.0
AIAM 3505C-H	3.5	5.0
AIAM 4501C-H	4.5	1.0
AIAM 4502C-H	4.5	2.0
AIAM 4503C-H	4.5	3.0
AIAM 4504C-H	4.5	4.0
AIAM 4505C-H	4.5	5.0



20 N.cm

**STRAIGHT CEMENTED  
ABUTMENT**

CODE	DIAM. (mm)	LENGTH (mm)
AIMP 3501C-H	3.5	1.0
AIMP 3502C-H	3.5	2.0
AIMP 3503C-H	3.5	3.0
AIMP 3504C-H	3.5	4.0
AIMP 3505C-H	3.5	5.0
AIMP 4501C-H	4.5	1.0
AIMP 4502C-H	4.5	2.0
AIMP 4503C-H	4.5	3.0
AIMP 4504C-H	4.5	4.0
AIMP 4505C-H	4.5	5.0

**TI BASE**

CODE	TRANSMUCOSAL LENGTH (mm)	CEMENTATION LENGTH (mm)
ICMT 0504	0.5	4.0
ICMT 0506	0.5	6.0
ICMT 2004	2.0	4.0
ICMT 2006	2.0	6.0
ICMT 3004	3.0	4.0
ICMT 3006	3.0	6.0

**LABORATORY SCREW****CODE**

PTMAML16	Hexagonal
PTL16	Titanium Hex.

**RETAINING SCREW****CODE**

PT 16	Titanium Hex.
	1.6mm screw

1.6mm screw

\*Hex Screw

\*Anti-Rotational Component

\*Squared Screw

\*Abutment Screw

# PROSTHETIC SEQUENCE 16°

## MULTI-UNIT ABUTMENT

MULTIPLE SCREW RETAINED RESTORATIONS



### IMPLANT

CODE EPIKUT S	CODE EPIKUT S PLUS	DIAM. (mm)	LENGTH (mm)
ILM 3585	ILM 3585N	3,5	8,5
ILM 3510	ILM 3510N	3,5	10
ILM 3511	ILM 3511N	3,5	11,5
ILM 3513	ILM 3513N	3,5	13
ILM 3515	ILM 3515N	3,5	15
ILM 3885	ILM 3885N	3,8	8,5
ILM 3810	ILM 3810N	3,8	10
ILM 3811	ILM 3811N	3,8	11,5
ILM 3813	ILM 3813N	3,8	13
ILM 3815	ILM 3815N	3,8	15
ILM 4085	ILM 4085N	4	8,5
ILM 4010	ILM 4010N	4	10
ILM 4011	ILM 4011N	4	11,5
ILM 4013	ILM 4013N	4	13
ILM 4015	ILM 4015N	4	15
ILM 4585	ILM 4585N	4,5	8,5
ILM 4510	ILM 4510N	4,5	10
ILM 4511	ILM 4511N	4,5	11,5
ILM 4513	ILM 4513N	4,5	13
ILM 4515	ILM 4515N	4,5	15
ILM 5085	ILM 5085N	5	8,5
ILM 5010	ILM 5010N	5	10
ILM 5011	ILM 5011N	5	11,5
ILM 5013	ILM 5013N	5	13
ILM 5015	ILM 5015N	5	15



20 N.cm

### STRAIGHT MULTI-UNIT ABUTMENT

CODE	DIAM. (mm)	LENGTH (mm)
MAM 4801 C	4.8	1
MAM 4802 C	4.8	2
MAM 4803 C	4.8	3
MAM 4804 C	4.8	4



20 N.cm

### INDEXED ANGLED MULTI-UNIT ABUTMENT

CODE	ANG.	DIAM. (mm)	LENGTH (mm)
MAAM 4802I	17°	4.8	2
MAAM 4803I	17°	4.8	3
MAAM 4804I	17°	4.8	4
MAAM 4832I	30°	4.8	2
MAAM 4833I	30°	4.8	3
MAAM 4834I	30°	4.8	4

\*Use hexagonal driver 1,2 mm



### OPEN TRAY TRANSFER

CODE
TMAM 4800

### ABUTMENT PROTECTOR

CODE
PMA 4855

5.0 mm profile



### CLOSED TRAY TRANSFER

CODE
TMFM 4800

**ANALOG**

**CODE**  
ANMA 4800

**TEMPORARY TITANIUM CYLINDER**

**CODE**

PTM 4800-3	For straight multi-unit
PTMS 4800-3	For straight multi-unit Suitable for laser welding
PTM 4800-2	For angled multi-unit

**CALCINABLE AND CO-CR CYLINDER**

**CODE**

CPM 4800-3	Plastic/For straight multi-unit
CLEM 4800-3	Cobalt chrome/For straight multi-unit
CPM 4800-2	Plastic/For angled multi-unit
CLEM 4800-2	Cobalt chrome/For angled multi-unit

**POLISHING PROTECTOR**

**CODE**  
PPM 01

**TI BASE**

**CODE**

PLAT (mm)	CEMENTATION LENGTH (mm)	
IMAT 02	4.8	2.0
IMAT 04	4.8	4.0
IMAT 06	4.8	6.0

**LABORATORY SCREW**

**CODE**

PL 1405 Short
PTMA 13-1 Long

**RETAINING SCREW**

**CODE** **HEIGHT (mm)**

PRH 20	2
--------	---

**RETAINING SCREW PACK 4**

**CODE** **HEIGHT (mm)**

PRH 30	3.0
--------	-----

**10 N.cm**

**10 N.cm**

**10 N.cm**

**10 N.cm**

- \*Hex Screw
- \*Anti-Rotational Component
- \*Squared Screw
- ◇ \*Abutment Screw




- Indicated for intraoral surgical placement in the maxilla, preferably in bones type III and IV (low density bones), for total edentulism cases, post extraction alveolus, immediate and delayed loading.
- It can be used in cases of total edentulous maxillae, especially in low density bones (bones type II and IV)
- The exclusive macro geometry guarantees precision and agility at the time of surgery.

#### INDICATIONS FOR CLINICAL USE:

- 3.8 - Anterior region
- 4.0 - Anterior and posterior region
- 4.5 - Posterior region

- Infra-bone installation
- Initial drill speed: 1200 rpm
- Speed of the drills 2.7 to 4.3mm: 800 rpm.
- Insertion speed: 20 to 40 rpm
- Maximum torque: 80 N.cm
- Immediate loading\*: recommended torque from 45 to 80 N.cm

\* Relative contraindication in patients with systemic or local problems and at the professional's discretion.



# EPIKUT S LONG DRILLING SEQUENCE

## FOR SOFT TYPE BONES

Drilling sequence  
used for bone type  
IV.

		1.200 RPM		800 RPM					
	∅ DIAM. (mm)	FL 2024 (A)	FHE 2324 (B)	FHE 3024 (C)	FHI 3324 (D)	FHI 3624 (E)	FHI 3824 (E+)	FHI 4024 (F)	FHI 4324 (G)
Epikut S Long	ILM38xx	3.8	●	●	●				
Epikut S Long Plus	ILM40xx	4.0	●	●	●	●			
Epikut S Long	ILM45xx	4.5	●	●	●	●	●		
Epikut S Long Plus									

## FOR MEDIUM TYPE BONES

Drilling sequence  
used for bone type  
II and III.

		1.200 RPM		800 RPM					
	∅ DIAM. (mm)	FL 2024 (A)	FHE 2324 (B)	FHE 3024 (C)	FHI 3324 (D)	FHI 3624 (E)	FHI 3824 (E+)	FHI 4024 (F)	FHI 4324 (G)
Epikut S Long	ILM38xx	3.8	●	●	●	●	●		
Epikut S Long Plus	ILM40xx	4.0	●	●	●	●	●	●	
Epikut S Long	ILM45xx	4.5	●	●	●	●	●	●	●
Epikut S Long Plus									

● USE OF DRILL IS OPTIONAL

## FOR HARD TYPE BONES

Drilling sequence  
used for bone type I.

1.200 RPM

800 RPM

	$\varnothing$ DIAM. (mm)	FL 2024 (A)	FHE 2324 (B)	FHE 3024 (C)	FHI 3324 (D)	FHI 3624 (E)	FHI 3824 (E+)	FHI 4024 (F)	FHI 4324 (G)
ILM38xx	3.8	•	•	•	•	•	•		
ILM40xx	4.0	•	•	•	•	•	•		
Epikut S Long	ILM45xx	4.5	•	•	•	•	•	•	•
Epikut S Long Plus									

## Technical measures



ILM 38xx  
ILM 38xxN



ILM 40xx  
ILM 40xxN



ILM 45xx  
ILM 45xxN



# EPIKUT S LONG PROSTHETIC SEQUENCE

## MULTI-UNIT ABUTMENT

MULTIPLE SCREW RETAINED RESTORATIONS



20 N.cm

### STRAIGHT MULTI-UNIT ABUTMENT

CODE	DIAM. (mm)	LENGTH (mm)
MAM 4801 C	4.8	1
MAM 4802 C	4.8	2
MAM 4803 C	4.8	3
MAM 4804 C	4.8	4



20 N.cm

### INDEXED ANGLED MULTI-UNIT ABUTMENT

CODE	ANG.	DIAM. (mm)	LENGTH (mm)
MAAM 4802I	17°	4.8	2
MAAM 4803I	17°	4.8	3
MAAM 4804I	17°	4.8	4
MAAM 4832I	30°	4.8	2
MAAM 4833I	30°	4.8	3
MAAM 4834I	30°	4.8	4

### IMPLANT

CODE EPIKUT S LONG	CODE EPIKUT S LONG PLUS	DIAM. (mm)	LENGTH (mm)
ILM 3818	ILM 3818N	3,8	18
ILM 3820	ILM 3820N	3,8	20
ILM 3822	ILM 3822N	3,8	22
ILM 3824	ILM 3824N	3,8	24
ILM 4018	ILM 4018N	4,0	18
ILM 4020	ILM 4020N	4,0	20
ILM 4022	ILM 4022N	4,0	22
ILM 4024	ILM 4024N	4,0	24
ILM 4518	ILM 4518N	4,5	18
ILM 4520	ILM 4520N	4,5	20



### OPEN TRAY TRANSFER

CODE  
TMAM 4800



ABUTMENT PROTECTOR  
CODE  
PMA 4855  
5.0 mm profile



### CLOSED TRAY TRANSFER

CODE  
TMFM 4800

\*Use hexagonal driver 1.2 mm.



**ANALOG**  
CODE  
ANMA 4800



10 N.cm

**TEMPORARY TITANIUM CYLINDER****CODE**

PTM 4800-3	For straight multi-unit
PTMS 4800-3	For straight multi-unit Suitable for laser welding
PTM 4800-2	For angled multi-unit



10 N.cm

**CALCINABLE AND CO-CR CYLINDER****CODE**

CPM 4800-3	Plastic/For straight multi-unit
CLEM 4800-3	Cobalt chrome/For straight multi-unit
CPM 4800-2	Plastic/For angled multi-unit
CLEM 4800-2	Cobalt chrome/For angled multi-unit

**POLISHING PROTECTOR**  
CODE PPM 01**LABORATORY SCREW**

CODE	DIAM. (mm)
PL 1405 Short	1.4
PTMA 13-1 Long	1.4



10 N.cm

**RETAINING SCREW**

CODE	HEIGHT (mm)
PRH 20	2



10 N.cm

**RETAINING SCREW PACK 4**

CODE	HEIGHT (mm)
PRH 30	3.0

◆ \*Hex Screw


◎ \*Anti-Rotational Component


■ \*Squared Screw

○ \*Abutment Screw

# EPIKUT SURGICAL KIT

MAXIMUM FUNCTIONALITY AND SIMPLICITY FOR YOUR SURGERIES







PRODUCT CODE: KCSE 02  
ORGANIZING BOX CODE: COSE 02

\*Check product availability in your country.

# EPIKUT S LONG SURGICAL KIT

MAXIMUM FUNCTIONALITY AND SIMPLICITY FOR YOUR SURGERIES





CODE: KCSEL 02i

ORGANIZING BOX CODE: COSEL 02

\*Check product availability in your country.

# EPIKUT SAFE DRILL KIT

## MAKING YOUR SURGERIES MORE PRACTICAL AND PRECISE

Performance and efficiency:  
exclusive polyacetal limiters with  
perfect fit and high resistance, which  
guarantees greater durability of the kit.

Bone Drill Limiters available  
for each drill diameter.


Prevent injuries to  
noble structures  
like nerves, maxillary  
sinus and nasal cavity.


The Epikut Safe Drill Kit  
is only compatible with the  
Epikut Surgical Kit.

Easy to use: color coding system,  
which facilitates clinical use.

Removable tray to  
facilitate cleaning.

For the Morse Taper installation to occur as recommended (infra-bone)  
it is necessary to use a limiter 1.5 mm greater than the desired depth.






CODE: KESD 02  
ORGANIZING BOX CODE: COESD 02

# PROSTHETIC KIT

## FUNCTIONAL, PRACTICAL AND COMPACT





\*Check product availability in your country.



CODE: KTMEC 02 | ORGANIZING BOX CODE: COTMEC 02

## DIGITAL SCREWDRIVERS







ITEM	CODE	DESCRIPTION	LENGTH	INDICATION
	CDQ 1224	SQUARE DIGITAL SCREWDRIVER 24.0MM	LONG	Used to set the square-fit locking screws (PTQ 2008, PTQH 18 and PTQ 2006). 1.3mm tip
	CDQ 1237	SQUARE DIGITAL SCREWDRIVER 37.0MM	EXTRA LONG	Used to set the square-fit locking screws (PTQ 2008, PTQH 18 and PTQ 2006). 1.3mm tip
	CLH 1277	HEX. SCREWDRIVER 77.0MM	EXTRA LONG	Lab screwdriver. Used to set retaining screws (PTL 16, PT 2006, PT 2008, PRH 20 and PRH 30) and lab screws. 1.2mm hexagonal tip
	CLQ 1277	HEX. SCREWDRIVER 77.0MM	EXTRA LONG	Lab screwdriver. Used to set the square-fit retaining screws (PTQ 2008, PTQH 18 and PTQ 2006). 1.3mm tip
	CRC 16	PROVISIONAL CYLINDER REMOVAL SCREWDRIVER	SHORT	Used to remove 1.6mm Cone Morse Strong SW provisional cylinder
	CRC 18	PROVISIONAL CYLINDER REMOVAL SCREWDRIVER	SHORT	Used to remove the 1.8 mm Cone Morse 11,5° provisional cylinder
	CDH 1620	HEX DIGITAL SCREWDRIVER 16MM	SHORT	Used to install the Multifunctional Abutment. 1.6mm Hexagonal Tip
	CDH 1624	HEX DIGITAL SCREWDRIVER 16MM	MEDIUM	Used to install the Multifunctional Abutment. 1.6mm Hexagonal Tip
	CCH 1620	HEX RATCHET WRENCH 16MM	SHORT	Used for the installation and torque of the Multifunctional Abutment. 1.6mm Hexagonal Tip
	CCH 1624	HEX RATCHET WRENCH 16MM	MEDIUM	Used for the installation and torque of the Multifunctional Abutment. 1.6mm Hexagonal Tip

## BONE PROFILING MILLING CUTTERS

ITEM	CODE	DESCRIPTION	INDICATION
	PO 4150	Platform 4.1 mm – External Hex.	Opens bone profile to 5.0 mm
	PO 5055	Platform 5.0 mm – External Hex.	Opens bone profile to 5.5 mm

\*Check product availability in your country.

## COUNTER-ANGLE SCREWDRIVER

ITEM	CODE	DESCRIPTION	LENGTH	INDICATION
	CTA 1224	ABUTMENT TORQUE SCREWDRIVER 24.0MM	LONG	Used to set the mini-abutment and conical abutment screw
	CTH 0924	COUNTER-ANGLE HEXAGONAL TORQUE SCREWDRIVER 24.0MM	LONG	Used for installation of Externa Hex. Tryon implant cover, two-pieces straight universal abut and angled universal abut.
	CTH 1220	COUNTER-ANGLE HEXAGONAL TORQUE SCREWDRIVER 20.0MM	SHORT	Used to set the mounting piece, healing, transfer, retaining screws (PTL 16, PT 2006, PT 2008, PRH 20 and PRH 30) and lab screws. 1.2mm hexagonal tip
	CTH 1224	COUNTER-ANGLE HEXAGONAL TORQUE SCREWDRIVER 24.0MM	LONG	Used to set the mounting piece, healing, transfer, retaining screws (PTL 16, PT 2006, PT 2008, PRH 20 and PRH 30) and lab screws. 1.2mm hexagonal tip
	CTH 1230	COUNTER-ANGLE HEXAGONAL TORQUE SCREWDRIVER 30.0MM	EXTRA LONG	Used to set the mounting piece, healing, transfer, retaining screws (PTL 16, PT 2006, PT 2008, PRH 20 and PRH 30) and lab screws. 1.2mm hexagonal tip
	CTHA 1220	ANGULAR MINI-ABUTMENT COUNTER-ANGLE HEXAGONAL TORQUE SCREWDRIVER 20.0MM	SHORT	Used to set the angular mini-abutment screw 1.2mm hexagonal tip (except for the Unitite angular mini-abutment).
	CTHA 1224	ANGULAR MINI-ABUTMENT COUNTER-ANGLE HEXAGONAL TORQUE SCREWDRIVER 24.0MM	LONG	Used to set the angular mini-abutment screw 1.2mm hexagonal tip (except for the Unitite angular mini-abutment).
	CTQ 20	SQUARE TORQUE SCREWDRIVER 20.0MM	SHORT	Used counter-angle to set square-fit retaining screws (PTQ 2008, PTQH 18 and PTQ 2006). 1.3mm tip
	CTQ 24	SQUARE TORQUE SCREWDRIVER 24.0MM	LONG	Used counter-angle to set square-fit retaining screws (PTQ 2008, PTQH 18 and PTQ 2006). 1.3mm tip
	CTQ 30	SQUARE TORQUE SCREWDRIVER 30.0MM	EXTRA LONG	Used counter-angle to set square-fit retaining screws (PTQ 2008, PTQH 18 and PTQ 2006). 1.3mm tip

\*Check product availability in your country.

## HELICAL MILLING CUTTERS

ITEM	CODE	MEASUREMENTS	DESCRIPTION
	FH 2010	ø 2,0x 10,0 mm	
	FH2020	ø 2,0x 18,0 mm	<ul style="list-style-type: none"><li>&gt; Surgical-grade stainless steel</li><li>&gt; Thermal treatment</li></ul>
	FH3010	ø 3,0x 10,0 mm	<ul style="list-style-type: none"><li>&gt; Laser markings</li><li>&gt; Used as a sequence to make the alveolus</li></ul>
	FH3020	ø 3,0x 18,0 mm	

## TREPHINE MILLING CUTTERS

ITEM	CODE	MEASUREMENTS	DESCRIPTION
	FTR 02	ø 2,0 mm	
	FTR04	ø 4,2 mm	<ul style="list-style-type: none"><li>&gt; Surgical-grade stainless steel</li><li>&gt; Thermal treatment</li></ul>
	FTR 05	ø 5,1 mm	<ul style="list-style-type: none"><li>&gt; Laser markings</li><li>&gt; May be used to remove implants, remove bone, and bone biopsy</li></ul>
	FTR 06	ø 6,1 mm	<ul style="list-style-type: none"><li>&gt; Measures refer to the inner diameter of the part</li></ul>
	FTR 08	ø 8,0 mm	

\*Check product availability in your country.

# MORE EASILY AND SAFETY FOR YOUR CLINICAL PROCEDURES

S.I.N. Implant System packaging is practical, maintaining the products in their integrity, facilitating the handling and the identification.



- › **01** The package is easy to open and handle even with gloves on.



- › **02** Transparency of package for optimal visibility of the implant.



- › **03** Separate compartments in same package for implant and cover.

- › **04** Snap-on top opening system ensures sterilization of the implant.



- › **05** With a proper connector, capture the implant with the counter angle key and move it until it reaches the perfect fit.



- › **06** The only implant system that offers the cover screw in the same packaging. To capture it, remove the cover screw from the tube cap and fit it on the 1.2 mm hexagonal digital key.

The implant should not be captured with the ratchet wrench.

# GENERAL INSTRUCTIONS

Special care and clarification on surgical instruments.



## CLEANING KIT CASE

- Remove manually all surgical instruments from the kit. Remove the kit box parts (lid, tray and bottom).
- Prepare the enzymatic detergent, according to manufacturer's recommendation.
- Immerse the trays into the prepared detergent solution and keep in contact for at least 5 minutes, then using a soft bristle brush, scrub the parts to remove organic matter from the products.
- Remove trays from detergent solution and rinse with tap water for 1 minute, repeat the rinse for two more times, a total of three rinses of 1 minute each.
- Visual inspection of each part for cleaning process residue or organic waste from product use.
- If residue is detected in the product, repeat the cleaning process until the residue is completely removed.
- Dry with a soft, clean, dry cloth or disposable paper.



## STERILIZATION

- Reusable Product and provided non-sterile.
- It must be clean and sterilized in autoclave before use.
- Dry all instruments before the steam sterilization cycle.
- The product must be enclosed in a steam sterilizable wrap.
- Steam sterilize in cycles of 121°C at 1 ATM pressure for 30 minutes or of 134°C at 2 ATM pressure for 20 minutes. Drying time 30 minutes.
- Always accommodate the case in autoclave over a plane surface and away of device walls.
- Never stack objects or other cases.



## CLEANING SURGICAL INSTRUMENTS

- Disassemble the product (if applicable). For the torque wrench, disassembly it completely, remove all the internal organic matter using tap water and follow to the next step only after performing such procedures.
- Prepare the enzymatic detergent according to the manufacturer's recommendation.
- Immerse all parts of the product into the prepared detergent solution and keep in contact for at least 5 minutes, then using soft bristle brush, scrub the parts to remove organic matter from the products.
- Remove parts from detergent solution and rinse with tap water for 1 minute, repeat the rinse for two more times, a total of three rinses of 1 minute each.
- Visual inspection of each part for cleaning process residue or organic waste from product use.
- If residue is detected in the product, repeat the cleaning process until the residue is completely removed.
- Dry with a soft, clean, dry cloth or disposable paper.
- Follow to sterilization process.

## CLEANING RECOMMENDATION

- Use the proper PPEs (gloves, masks, goggles, caps, etc.).
- Start the cleaning right after the surgical use.
- Never let the instruments dry with organic waste after the surgical use.
- Never let the instrument dry naturally after cleaning.
- Never use saline solutions, include sodium hypochlorite, disinfectant, hydrogen peroxide or alcohol for cleaning or rinsing the surgical instruments and Kits.
- Never use steel wool and abrasive products, so that the instruments are not damaged.
- Do not stack the instruments in lots to avoid the deformation of smaller and delicate pieces.

## STERILIZATION RECOMMENDATIONS

- Sterilize the products in the same day or one day earlier the procedure.
- The chemical sterilization is not recommended, once some products may cause the discoloration and damages to the case.
- Do not use temperature higher than 60°C to drying process.
- Do not use dry heat stoves for sterilization of the instruments and kits from S.I.N.

Discover **IMPLANTAT**,  
the educational habitat of S.I.N. Implant System.  
An online teaching platform created to make more professionals  
accelerate their career and increase their success.

Access  
**IMPLANTAT.GLOBAL**  
or scan the QR Code and begin  
your journey of knowledge now!





Discover **Implantat**, the educational habitat of S.I.N. Implant System.

 [implantat.global](https://implantat.global)

1-866-507-9315

[SINdentalusa.com](http://SINdentalusa.com) | [sinimplantsystem.com](http://sinimplantsystem.com)

---

Visit our social networks: