



• SMOOTH AND POLISHED SURFACE IN THE CERVICAL AREA:

Coated with a smooth surface, facilitating easy cleaning, and reducing biofilm buildup.

• AGGRESSIVE AND COMPRESSIVE THREADS:

Offers unique performance and enables faster implant locking.

• 16° MORSE TAPER CONNECTION:

Generates high mechanical strength.

• SHARP APEX TO FACILITATE IMPLANT INSERTION:

Provides optimal results for severe cases of maxillary resorption.



# IDEAL FOR IMMEDIATE LOADING

#### LIST OF APPLICABLE PRODUCTS CODE THREAD DIAMETER (MM) LENGTH (MM) PLATFORM DIAMETER (MM) ILMZ 4030EN Ø4.0 Ø4.0 30 ILMZ 4032EN Ø4.0 040 32.5 ILMZ 4035EN Ø4.0 35 Ø4.0 Ø4.0 ILMZ 4037EN Ø4.0 37.5 II M7 4040FN 040 040 40 42.5 ILMZ 4042EN 040 Ø4.0 ILMZ 4045EN Ø4.0 Ø4.0 45 II M7 4047FN Ø4 0 Ø4 0 47.5 50 ILMZ 4050EN Ø4.0 Ø4.0 II M7 4052FN Ø4 0 Ø4 0 525 ILMZ 4055EN 040 040 55 II M7 4057FN Ø4 0 Ø4 0 57.5

Ø4.0

60

Ø4.0

ILMZ 4060EN

## Immediate results and shorter treatment time

The Zygomatic implant offers an immediate solution to patients as an alternative to bone graft procedures. It ensures stability, predictability, and treatment precision, meeting patient expectations through its versatile design and options tailored to anatomical features and tissue relationships.




# PROSTHETIC COMPONENTS


Mult-unit conical abutment types with angulations of 17°, 30°, 45°, 45° Slim,  $52^{\circ}$ , and  $60^{\circ}$  with gingival heights ranging from 2.0 to 4.0 mm. Made of titanium alloy and anodized in gold color, with a locking torque of 20Ncm.


Protectors are screwed over the standard multi-unit abutments with manual torque (maximum 10ncm) and should remain in position until the provisional or definitive prosthesis is manufactured.











CLEM 4800-2 Cobalt chrome for angled multi-unit



### Discover the S.I.N. 360 Advanced Digital Workflow



7755 E Gray Rd Scottsdale, AZ 85260

(866) 507-9315

www.sin360.us

Visit our social channels:

Instagram
@sindentalusa

Facebook
@sindentalusa

**YouTube**@sindentalusa