

Versalis


 Versalis

Condensed Catalog: Rev 01.0

S.I.N. 360 

A SMART COMBINATION OF DESIGN AND VERSATILITY.

PLUS (HAnano)



"PLUS" SURFACE




› Hybrid Design, Smarter Mechanics:

A dual-tapered body that balances compression and stress distribution, offering better bone-implant contact and long-term success, enhanced by the ultra-thin Surface Plus, treated with double acid-etching and HAnano coating.



› An Implant with Diverse Possibilities:

Morse Taper 16° connection making your clinical day-to-day easier.



› Exclusive Plus Surface:


Developed in the main universities of Sweden, the Plus surface which is produced by double acid-etching followed by application of a hydroxyapatite coating HAnano, proven by over 50 preclinical studies.



› From 8.5 to 24 mm. One Surgical Kit:

Standard or long, all lengths installed with the same kit.

Double Acid Etch



"DAE" SURFACE




› Easy Bicortical Engagement

Dual taper body and aggressive apical design ensure easier bicortical anchorage, even in challenging bone anatomy.



› Dual Thread:

Faster insertion, stable from the start. Less time placing, more time producing.




› Manufactured in Cold Worked Grade IV Titanium:

Super light metal, very resistant to corrosion, wear and fracture.



› Cutting Apex for Maximum Primary Stability:


No extra steps. No over-prep. Versalis cuts its own path and achieves high initial stability with minimal bone trauma.



› Full Surface Treatment:


Double acid etching on the entire surface.

ONE CONNECTION PLATFORM



All Versalis Implants feature the 16° Morse Taper Connection with Platform Switch, providing one single connection for all implant diameters.

TECHNICAL MEASURES



DRILLING SEQUENCE: VERSALIS LONG

BONE TYPES III AND IV




Ø3.5	●	●	● ● ●		
Ø3.8	●	●	●	● ●	
Ø4.0	●	●	●	●	● ●


● Optional ● Optional in case of very soft bone Initial Drill speed of 1200 rpm Conical Drill speed of 800 rpm

SINGLE CONNECTION SMARTER WORKFLOW

All Versalis Long Implants feature the 16° Morse Taper Connection with Platform Switch, providing one single connection for all implant diameters



TECHNICAL MEASURES




MT 16° PROSTHETIC SEQUENCE

DIRECT SEQUENCE OVER THE IMPLANT (DIGITAL)


SINGLE RESTORATIONS




IMPLANT			
CODE VERSALIS S	CODE VERSALIS S PLUS	DIAM. (mm)	LENGTH (mm)
IVMT 3585	IVMT 3585N	3.5	8.5
IVMT 3510	IVMT 3510N	3.5	10.0
IVMT 3511	IVMT 3511N	3.5	11.5
IVMT 3513	IVMT 3513N	3.5	13.0
IVMT 3515	IVMT 3515N	3.5	15.0
IVMT 3518	IVMT 3518N	3.5	18.0
IVMT 3885	IVMT 3885N	3.8	8.5
IVMT 3810	IVMT 3810N	3.8	10.0
IVMT 3811	IVMT 3811N	3.8	11.5
IVMT 3813	IVMT 3813N	3.8	13.0
IVMT 3815	IVMT 3815N	3.8	15.0
IVMT 3818	IVMT 3818N	3.8	18.0
IVMT 4085	IVMT 4085N	4.0	8.5
IVMT 4010	IVMT 4010N	4.0	10.0
IVMT 4011	IVMT 4011N	4.0	11.5
IVMT 4013	IVMT 4013N	4.0	13.0
IVMT 4015	IVMT 4015N	4.0	15.0
IVMT 4018	IVMT 4018N	4.0	18.0
IVMT 4385	IVMT 4385N	4.3	8.5
IVMT 4310	IVMT 4310N	4.3	10.0
IVMT 4311	IVMT 4311N	4.3	11.5
IVMT 4313	IVMT 4313N	4.3	13.0
IVMT 4315	IVMT 4315N	4.3	15.0
IVMT 4318	IVMT 4318N	4.3	18.0
IVMT 5085	IVMT 5085N	5.0	8.5
IVMT 5010	IVMT 5010N	5.0	10.0
IVMT 5011	IVMT 5011N	5.0	11.5
IVMT 5013	IVMT 5013N	5.0	13.0
IVMT 5015	IVMT 5015N	5.0	15.0
IVMT 5018	IVMT 5018N	5.0	18.0
IVMT 6085	IVMT 6085N	6.0	8.5
IVMT 6010	IVMT 6010N	6.0	10.0
IVMT 6011	IVMT 6011N	6.0	11.5
IVMT 6013	IVMT 6013N	6.0	13.0
IVMT 7085	IVMT 7085N	7.0	8.5
IVMT 7010	IVMT 7010N	7.0	10.0
IVMT 7011	IVMT 7011N	7.0	11.5
IVMT 7013	IVMT 7013N	7.0	13.0



TITANIUM HEALING CAP		
CODE	DIAM. (mm)	HEIGHT (mm)
CIM 3502C	3.5	2.0
CIM 3504C	3.5	4.0
CIM 3506C	3.5	6.0
CIM 4502C	4.5	2.0
CIM 4504C	4.5	4.0
CIM 4506C	4.5	6.0



PEEK HEALING CAP		
CODE	PROFILE DIAM. (mm)	HEIGHT (mm)
CPCM 0504	5.0	4.0
CPCM 0804	8.0	4.0
CPCM 0508	5.0	8.0
CPCM 0808	8.0	8.0



MT 16° SCANNING JIG	
CODE	
JBSWCM	



MT 16° SCANNING JIG	
CODE	
JBSWCMC	



DIGITAL ANALOG - MT 16°	
CODE	
ADCM	



TITANIUM INTERFACE MT 16° SIRONA				
S.I.N. PLATFORM	DESCRIPTION	TRANSMUCOSAL HEIGHT (mm)	LENGTH (mm)	SIRONA LIBRARY
ICM 0804	0.8x4.0	0.8	4.0	ATOS 3.5/4.0 – ATOS 4.5/5.0
ICM 2004	2.0x4.0	2.0	4.0	ATOS 3.5/4.0 – ATOS 4.5/5.0



TITANIUM INTERFACE MT 16°			
CODE	DESCRIPTION	TRANSMUCOSAL HEIGHT (mm)	LENGTH (mm)
ICMT 0504	0.5X4	0.5	4.0
ICMT 0506	0.5X6	0.5	6.0
ICMT 2004	2.0X4	2.0	4.0
ICMT 2006	2.0X6	2.0	6.0
ICMT 3004	3.0X4	3.0	4.0
ICMT 3006	3.0X6	3.0	6.0



CUSTOM PREMILL ABUTMENTS	
CODE	DESCRIPTION
AFTCM 1601DS	16° MT Premill Narrow (DESS)
AFTCM 1602DS	16° MT Premill Wide (DESS)
AFTCM 1601 - MK	11.5° Premill (Medentika)
AFTCM 1602 - MK	15.8° Premill (Medentika)

DRIVERS

Driver Handpiece
Hex. 1.2mm Short
(CTH 1220)

Driver Handpiece
Hex. 1.2mm Medium
(CTH 1224)

Driver Handpiece
Hex. 1.2mm Long
(CTH 1230)

Driver Ratchet
Hex. 1.2mm Short
(CDHC 20)

Driver Ratchet
Hex. 1.2mm
(CDHC 24)


*Check product availability in your country.

- * Analog sequence
- * Digital sequence
- * Hex driver
- * Anti-Rotational component
- * Squared Screw
- * Abutment Screw
- * Rotational component

MT 16° PROSTHETIC SEQUENCE

DIRECT SEQUENCE OVER THE IMPLANT (ANALOG)

SINGLE RESTORATIONS




*Check product availability in your country.

MT 16° PROSTHETIC SEQUENCE

MULTI-UNIT ABUTMENT (ANALOG & DIGITAL)

MULTIPLE SCREW RETAINED RESTORATIONS



*Check product availability in your country.

VERSALIS SURGICAL KIT ORGANIZING BOX



1 LANCE DRILL Ø2.0MM (FRL 20)	9 CONICAL DRILL Ø7.0MM SOFT (FRCV 70S)	17 COUNTERSINK DRILL Ø3.5MM (FRCV 35)	25 DIRECTION INDICATOR 3.8mm (ID 38)	33 LANCE DRILL Ø2.0MM LONG (FRL 20L)	41 DIGITAL DRIVER HEX.1.2 MEDIUM (CDH 1224)
2 HELICAL DRILL Ø2.0MM (FRH 20)	10 LINDEMANN DRILL (FRL 01)	18 COUNTERSINK DRILL Ø3.8MM (FRCV 38)	26 DIRECTION INDICATOR 4.0mm (ID 40)	34 CONICAL DRILL Ø2.0MM LONG (FRL 20L)	42 HANDPIECE DRIVER ABUT. Ø1.2MM MEDIUM (CTA 1224)
3 CONICAL DRILL Ø3.5MM SOFT (FRCV 35S)	11 CONICAL DRILL Ø3.5MM HARD (FRCV 35H)	19 COUNTERSINK DRILL Ø4.0MM (FRCV 40)	27 DIRECTION INDICATOR 4.3mm (ID 43)	35 ANGLED DIRECTION INDICATOR 17° CONICAL DRILL Ø3.8MM LONG (FRCVL 38)	43 SURGICAL TORQUE RATCHET (TMECC 01EN)
4 CONICAL DRILL Ø3.8MM SOFT (FRCV 38S)	12 CONICAL DRILL Ø3.8MM HARD (FRCV 38H)	20 COUNTERSINK DRILL Ø4.3MM (FRCV 43)	28 DIRECTION INDICATOR 5.1mm (ID 50)	36 CONICAL DRILL Ø3.8MM LONG (FRCVL 38)	44 SURGICAL TORQUE RATCHET (TMECC 02)
5 CONICAL DRILL Ø4.0MM SOFT (FRCV 40S)	13 CONICAL DRILL Ø4.0MM HARD (FRCV 40H)	21 COUNTERSINK DRILL Ø5.0MM (FRCV 50)	29 HANDPIECE DRIVER IMP. MT16° SHORT (IDMT 01C)	37 CONICAL DRILL Ø4.0MM LONG (FRCVL 40)	45 DEPTH PROBE STANDARD/LONG (SOP 02)
6 CONICAL DRILL Ø4.3MM SOFT (FRCV 43S)	14 CONICAL DRILL Ø4.3MM HARD (FRCV 43H)	22 LANCE DRILL Ø2.0MM LONG (FRL 20L)	30 HANDPIECE DRIVER IMP. MT16° SHORT (IDMT 01C)	38 CONICAL DRILL Ø4.3MM LONG (FRCVL 43)	
7 CONICAL DRILL Ø5.0MM SOFT (FRCV 50S)	15 CONICAL DRILL Ø5.0MM HARD (FRCV 50H)	23 HANDPIECE DRIVER IMP. MT16° LONG (IDMT 01L)	31 ANGLED DIRECTION INDICATOR 30° (IDMT 01C)	39 DIRECTION INDICATOR OF DRILL (IDMT 03)	
8 CONICAL DRILL Ø6.0MM SOFT (FRCV 60S)	16 RATCHET DRIVER IMP. MT SHORT (IDMT 01C)	24 DIRECTION INDICATOR 3.5mm (ID 35)	32 HANDPIECE DRIVER IMP. MT16° LONG (IDMT 01L)	40 HANDPIECE DRIVER ABUT. Ø1.2MM MEDIUM (CTH 1224)	

PROSTHETIC KIT ORGANIZING BOX



1 Torque Wrench (TMEC)	6 Squared Driver 1.3x20mm (CQTM 20)
2 Hex. Driver 0.9x20mm (CCH 0920)	7 Squared Driver 1.3x24MM (CQTM 24)
3 Hex. Driver 1.2x24mm (CDHC 24)	8 Multi-Unit/Conical Driver (CDAC 20)
4 Hex. Driver 1.2x20mm (CDHC 20)	9 Angled Multi-Unit Driver 1.2mm (CHTMA 24)
5 Hex. Driver 1.6mm (CCH 1624)	10 O-Ring Driver (CCAO 20)
	11 Digital Adapter (CPQ 02)

Remote Anchorage



Epikut Long RA

Our **Epikut[®] Long RA** is a specialized dental implant designed for predicatable stability. Its smooth, machined collar combines dual thread technology and a narrow apex to make this one of our most advanced implants.

Part Numbers

DAE: ILMT xxxx

Length (mm)	Ø Diameter (mm)			
	3.8	4.0	4.2	4.5
18	ILMT 3818	ILMT 4018	ILMT 4218	ILMT 4518
20	ILMT 3820	ILMT 4020	ILMT 4220	ILMT 4520
22	ILMT 3822	ILMT 4022	ILMT 4222	ILMT 4522
24	ILMT 3824	ILMT 4024	ILMT 4224	ILMT 4524



KCSEP 01

Remote Anchorage

Zygomatic


Our **Zygomatic** implants are designed with a smooth polished surface in the cervical area, aggressive and compressive threads, with a sharp apex to facilitate implant insertion; providing superior clinical results for severe cases of maxillary resorption.

Part Numbers

DAE: ILMZ xxxxE
Plus (HAnano): ILMZ xxxxEEN

Diameter: Ø4.0 mm

Length (mm)	Part Number (ILMZ xxxxE)	Length (mm)	Part Number (ILMZ xxxxE)	Length (mm)	Part Number (ILMZ xxxxE)
30	ILMZ 4030E	40	ILMZ 4040E	50	ILMZ 4050E
32.5	ILMZ 4032E	42.5	ILMZ 4042E	52.5	ILMZ 4052E
35	ILMZ 4035E	45	ILMZ 4045E	55	ILMZ 4055E
37.5	ILMZ 4037E	47.5	ILMZ 4047E	57.5	ILMZ 4057E
				60	ILMZ 4060E



KZ 01



Shop now

S.I.N. 360 

1-866-507-9315

sin360.us | shop.sin360.us



Visit our social networks



@sindentalUSA